


Foods

Flour quality control

With a TAGARNO microscope



Optimize workflows

Counting specs in no time

With all settings saved as a template, tedious routines have been reduced to one click

The TAGARNO Speck counter app optimizes this workflow to the extent that the delighted operator will find it all flow and almost no work, while at the same time delivering much improved quality in the form of objective and reproducible analyses.

The operator simply places the sample under the microscope and clicks to run the analysis. They usually do not need adjust the settings – such as choosing the minimum size of the speck or adjusting the lightness filter – because they did this the first time, and then saved all the settings as a template that can easily be recalled. Instead, one click is all it takes to count every single speck.

Yesteryear's manual counting and noting down is now fully automated and instantaneous. The result is not only a much more efficient process, but also perfectly consistent inspections.

“**Counting speck and bran
in a flour sample is easier
than ever before**”

The perfect package for

Flour quality control

Microscope



Accessories

Lens +10
Magnify at the right level

Ring light white
View natural sample colors

Foot switch
Control zoom levels with feet

Apps



Speck counter app
Count sample impurities

“Counting speck and bran
in a flour sample is **easier**
than ever before”

Contact

**Want to
know more?**

tagarno.com

TAGARNO