



Image comparison app

Build image library with own photos

Compare products to your own reference photos. Choose side-by-side view, alternate display or split view mode.

“With reference photos, we can easily spot **deviations** in seed samples. It’s great

TAGARNO

Results

Increase efficiency by comparing

Increase operator efficiency by using photos of seed samples as a visual guide instead of relying on memory. Add new photos to the library at any point.

How to

Capture, save and recall photos



Capture reference photo

- ✓ Save a photo of current sample
- ✓ Move the reference indicators to a suitable location
- ✓ If needed, crop the photo to only save relevant area
- ✓ Save the photo as a template for later use



Recall photos

- ✓ Select a template/reference photo from the dropdown menu
- ✓ Choose how you want the reference photo displayed
- ✓ Delete template when it's no longer needed
- ✓ [Read more about Image comparison app](#)

TAGARNO

tagarno.com