



Measurement app

# Measure directly on the microscope

No more rulers. Measure length, diameter and more directly on the microscope. Great for tracking seed growth.

“With onscreen measurements, we can easily track and document seed growth”

**TAGARNO**

## Results

# Increase accuracy and efficiency

Choose one of the many measurement options to track seed growth. Afterwards, save a photo of your work for internal documentation or share it externally.

## How to

# Measure, save and share



## Measure

- ✓ Use tool menu to choose your preferred measurement
- ✓ Choose between length, polyline, area, diameter and more
- ✓ Add text boxes and arrows to highlight an area of the inspection photo



## Save and share

- ✓ Save an image of your inspection with or without the measurements you've made
- ✓ Send the photo to colleagues or supply chain partners
- ✓ [Read more about Measurement app](#)

**TAGARNO**

[tagarno.com](https://tagarno.com)