



Ruler app

Measure objects for easy checks

The integrated Ruler app, a piece of measurement software for digital microscopes, makes sure that measuring any magnified object is quick and precise.



On-screen measurement makes it
easy to do a **quick quality check**

TAGARNO

Optimize

Activate ruler on screen and go

Simply activate the Ruler graphics vertically and/or horizontally as an overlay directly on the screen. Then place the object and change the color of the ruler to create a contrast to the object color. No external computer needed, the microscope is the computer. Ready to inspect and quality assure, just like that.

How to

Easy and colorful



Easy but custom



Effortless measurements

Painless and accurate measurements of any object.



Choose direction

Take a pick between having a ruler displayed horizontally or vertically. It's also possible to trigger both rulers at the same time.



Choose the color



Change color

Simply change the color of the ruler to create a contrast to the object.



Precision

Precision of ruler: 0.5 %.

TAGARNO

tagarno.com