



DIGITAL MICROSCOPE

TAGARNO ZIP

Compact and light weight

“We can perform twice as many
analyses as before. **We save time**”

Jens Asmussen, Laboratory Manager
DLG, Denmark

TAGARNO



TAGARNO ZIP

Effective quality control

Using a digital microscope makes everything stand out.

With a digital microscope, you can view your seed and grain sample test in extremely high image quality, no matter the magnification level.

The clear images make it easier to perform processes such as seed germination testing, variety identification and purity determination. Foreign varieties, mites or fungal diseases in your sample will also be easy to spot.

The setup with a monitor is also more comfortable and less tiring for the operator, creating healthier and happier employees.

Lastly, by connecting the microscope to a computer and taking a screenshot of your work, you will have image documentation readily at hand.

“Performing visual identification processes has become much **more efficient**”

Helle Foldgaard, laboratory technician
Nordic Seed, Denmark

TAGARNO ZIP

A simple digital microscope

Despite its sleek and simple appearance, ZIP produces ultra sharp images at all magnification levels and with no delay. This creates excellent conditions for any visual inspection task.



Rotate the camera

For extra flexibility, the camera head can rotate more than 180 degrees.

Unfold and you're off

With only one height setting, there's no fiddling with focus heights here.

Lots of room

There's plenty of room beneath the camera to do your inspection.

Small footprint

Being only 155 mm wide, TAGARNO ZIP doesn't take up much space.

Tech specs

Take a closer look

The microscope might be small, but don't underestimate the power of the ZIP.

Magnification levels	1.7 – 53x	+
Image quality	FHD/1080p	+
Frame rates	60FPS	+
Optical zoom	30x	+
Working depth	280 mm	+
Working height	245 mm	+
Field of view	11 – 290 mm	+

53x

Magnification level

60 FPS

Ultra-high frame rates

280MM

Working depth

1080P

Image quality

The perfect package

Everything you need

Mix and match. Add any accessory to your bundle and you are good to go.

Microscope



Accessories

XKEY control box
Use microscope features

Lens +4
Magnify at the right level

Glass table
Remove shadows

Contact

Want to know **more?**

Take a look at tagarno.com or contact your local TAGARNO distributor if you want to see how you can optimize processes and improve quality.

tagarno.com

TAGARNO